

Supplementary Material

An efficient debromination technique using PMHS with a number of ligands containing different functional groups

Md Yeunus Mian,^a Prithu Mondal,^a Dishary Sharmin,^a Kamal P. Pandey,^a Farjana Rashid,^a Sepideh Rezvanian,^a Lalit K. Golani,^a V.V.N. Phani Babu Tiruveedhula,^a Vilashini Rajaratnam¹, Shama P. Mirza¹, John D. Chan², Jeffrey M. Witkin¹, and James M. Cook*¹

¹Department of Chemistry and Biochemistry, Milwaukee Institute of Drug Discovery, University of Wisconsin-Milwaukee, Milwaukee, WI-53211



²Department of Chemistry, University of Wisconsin-Oshkosh, Oshkosh, WI-54901


E-mail: capncook@uwm.edu

Table of Contents

| | |
|--|-----|
| 1. ¹ HNMR spectra of compound 2 | S3 |
| 2. ¹³ CNMR spectra of compound 2 | S3 |
| 3. High resolution mass spectra of compound 2 | S4 |
| 4. ¹ HNMR spectra of compound 3 | S5 |
| 5. ¹³ CNMR spectra of compound 3 | S5 |
| 6. High resolution mass spectra of compound 3 | S6 |
| 7. ¹ HNMR spectra of compound 5b | S7 |
| 8. ¹³ CNMR spectra of compound 5b | S7 |
| 9. High resolution mass spectra of compound 5b | S8 |
| 10. ¹ HNMR spectra of compound 5c | S9 |
| 11. ¹³ CNMR spectra of compound 5c | S9 |
| 12. High resolution mass spectra of compound 5c | S10 |
| 13. ¹ HNMR spectra of compound 5d | S11 |
| 14. ¹³ CNMR spectra of compound 5d | S11 |
| 15. High resolution mass spectra of compound 5d | S12 |
| 16. ¹ HNMR spectra of compound 5e | S13 |
| 17. ¹³ CNMR spectra of compound 5e | S13 |
| 18. High resolution mass spectra of compound 5e | S14 |
| 19. ¹ HNMR spectra of compound 5f | S15 |
| 20. ¹³ CNMR spectra of compound 5f | S15 |
| 21. High resolution mass spectra of compound 5f | S16 |
| 22. ¹ HNMR spectra of compound 5g | S17 |

| | |
|--|-----|
| 23. $^{13}\text{CNMR}$ spectra of compound 5g | S17 |
| 24. High resolution mass spectra of compound 5g | S18 |
| 25. $^1\text{H}\text{NMR}$ spectra of compound 5h | S19 |
| 26. $^{13}\text{CNMR}$ spectra of compound 5h | S19 |
| 27. High resolution mass spectra of compound 5h | S20 |
| 28. $^1\text{H}\text{NMR}$ spectra of compound 7a | S21 |
| 29. $^{13}\text{CNMR}$ spectra of compound 7a | S21 |
| 30. High resolution mass spectra of compound 7a | S22 |
| 31. $^1\text{H}\text{NMR}$ spectra of compound 7b | S23 |
| 32. $^{13}\text{CNMR}$ spectra of compound 7b | S23 |
| 33. High resolution mass spectra of compound 7b | S24 |
| 34. $^1\text{H}\text{NMR}$ spectra of compound 7b | S25 |
| 35. $^{13}\text{CNMR}$ spectra of compound 7b | S25 |
| 36. High resolution mass spectra of compound 7b | S26 |
| 37. $^1\text{H}\text{NMR}$ spectra of compound 7d | S27 |
| 38. $^{13}\text{CNMR}$ spectra of compound 7d | S27 |
| 39. High resolution mass spectra of compound 7d | S28 |
| 40. $^1\text{H}\text{NMR}$ spectra of compound 9a | S29 |
| 41. High resolution mass spectra of compound 9a | S30 |
| 42. Debromination on a double bond containing substrate | S31 |

5-phenyl-1,3-dihydro-2H-benzo[e][1,4]diazepin-2-one(2)



Formula Predictor Report • COMP 2_12082020 Analysis_49.lcd

Page 1 of 1

Data File: C:\LabSolutions\Dat\ Vilashini Rajaratnam\12082020 Analysis\COMP 2_12082020 Analysis_49.lcd

| Elmt | Val. | Min | Max | Use Adduct |
|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------------|
| H | 1 | 10 | 15 | N | 3 | 1 | 3 | Si | 4 | 0 | 0 | Br | 1 | 0 | 0 | H |
| 2H | 1 | 0 | 0 | O | 2 | 1 | 3 | S | 2 | 0 | 0 | I | 3 | 0 | 0 | K |
| C | 4 | 10 | 25 | F | 1 | 0 | 1 | Cl | 1 | 0 | 0 | | | | | |

Error Margin (ppm): 300

DBE Range: -100.0 - 2000.0

Electron Ions: both

HC Ratio: unlimited

Apply N Rule: no

Use MSn Info: yes

Max Isotopes: all

Isotope RI (%): 1.00


Isotope Res: 10000

MSn Iso RI (%): 75.00


MSn Logic Mode: AND

Max Results: 10


Event#: 1 MS(E+) Ret. Time : 0.893 > 0.907 > 1.147 > 1.154 Scan# : 135 > 137 > 173 > 175






Measured region for 237.1007 m/z



C15 H12 N2 O [M+H]+ : Predicted region for 237.1022 m/z



| Rank | Score | Formula (M) | Ion | Meas. m/z | Pred. m/z | Df. (mDa) | Df. (ppm) | Iso | DBE |
|------|-------|--------------|--------|-----------|-----------|-----------|-----------|-------|------|
| 1 | 27.96 | C15 H12 N2 O | [M+H]+ | 237.1007 | 237.1022 | +1.5 | +6.33 | 36.45 | 11.0 |

5-(2-chlorophenyl)-1,3-dihydro-2H-benzo[e][1,4]diazepin-2-one(3)

Formula Predictor Report - COMP 8_12082020 Analysis_37.lcd

Page 1 of 1

Data File: C:\LabSolutions\Data\Vilashini Rajaratnam\12082020 Analysis\COMP 8_12082020 Analysis_37.lcd

| Elmt | Val. | Min | Max | Use Adduct |
|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------------|
| H | 1 | 10 | 30 | N | 3 | 1 | 3 | Si | 4 | 0 | 0 | Br | 1 | 0 | 0 | H |
| 2H | 1 | 0 | 0 | O | 2 | 1 | 3 | S | 2 | 0 | 0 | I | 3 | 0 | 0 | K |
| C | 4 | 10 | 25 | F | 1 | 0 | 1 | Cl | 1 | 0 | 1 | | | | | |

Error Margin (ppm): 300

HC Ratio: unlimited

Max Isotopes: all

MSn Iso RI (%): 75.00

DBE Range: +100.0 - 2000.0

Apply N Rule: no

Isotope RI (%): 1.00

MSn Logic Mode: AND


Electron Ions: both

Use MSn Info: yes


Isotope Res: 10000

Max Results: 10


Event#: 1 MS(E+) Ret. Time : 0.920 > 1.053 > 0.467 > 0.598 Scan# : 139 > 159 > 71 > 91






Measured region for 271.0608 m/z



C15 H11 N2 O Cl [M+H]+ : Predicted region for 271.0633 m/z



| Rank | Score | Formula (M) | Ion | Meas. m/z | Pred. m/z | Df. (mDa) | Df. (ppm) | Iso | DBE |
|------|-------|-----------------|--------|-----------|-----------|-----------|-----------|-------|------|
| 3 | 31.61 | C15 H11 N2 O Cl | [M+H]+ | 271.0608 | 271.0633 | -2.5 | -9.22 | 66.12 | 11.0 |

(R)-7-bromo-3-methyl-5-phenyl-1,3-dihydro-2H-benzo[e][1,4]diazepin-2-one(5b)

Formula Predictor Report - COMP 4_12082020 Analysis_45.lcd

Page 1 of 1

Data File: C:\LabSolutions\Dat\Viashini Rajaratnam\12082020 Analysis\COMP 4_12082020 Analysis_45.lcd


| Elmt | Val. | Min | Max | Use Adduct |
|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------------|
| H | 1 | 10 | 30 | N | 3 | 1 | 3 | Si | 4 | 0 | 0 | Br | 1 | 0 | 0 | H |
| 2H | 1 | 0 | 0 | O | 2 | 1 | 3 | S | 2 | 0 | 0 | I | 3 | 0 | 0 | K |
| C | 4 | 10 | 25 | F | 1 | 0 | 0 | Cl | 1 | 0 | 0 | | | | | |

Error Margin (ppm): 300
 HC Ratio: unlimited
 Max Isotopes: all
 MSn Iso RI (%): 75.00


DBE Range: -100.0 -> 2000.0
 Apply N Rule: no
 Isotope RI (%): 1.00
 MSn Logic Mode: AND

Electron Ions: both
 Use MSn Info: yes
 Isotope Res: 10000
 Max Results: 10


Event#: 1 MS(E+) Ret. Time : 1.480 -> 1.493 -> 0.693 -> 0.760 Scan# : 223 -> 225 -> 105 -> 115



Measured region for 251.1155 m/z






C16 H14 N2 O [M+H]+ : Predicted region for 251.1179 m/z



| Rank | Score | Formula (M) | Ion | Meas. m/z | Pred. m/z | Df. (mDa) | Df. (ppm) | Iso | DBE |
|------|-------|--------------|--------|-----------|-----------|-----------|-----------|-------|------|
| 1 | 25.99 | C16 H14 N2 O | [M+H]+ | 251.1155 | 251.1179 | +2.4 | +9.56 | 58.54 | 11.0 |

(S)-7-bromo-3-methyl-5-phenyl-1,3-dihydro-2H-benzo[e][1,4]diazepin-2-one(5c)



Formula Predictor Report - COMP 9_12082020 Analysis_35.lcd

Page 1 of 1

Data File: C:\LabSolutions\LabSolutions\12082020 Analysis\COMP 9_12082020 Analysis_35.lcd

| Elmt | Val. | Min | Max | Use Adduct |
|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------------|
| H | 1 | 10 | 30 | N | 3 | 1 | 3 | Si | 4 | 0 | 0 | Br | 1 | 0 | 0 | H |
| 2H | 1 | 0 | 0 | O | 2 | 1 | 3 | S | 2 | 0 | 0 | I | 3 | 0 | 0 | K |
| C | 4 | 10 | 25 | F | 1 | 0 | 1 | Cl | 1 | 0 | 1 | | | | | |

Error Margin (ppm): 300

HC Ratio: unlimited

Max Isotopes: all

MSn Iso RI (%): 75.00

DBE Range: 100.0 - 2000.0

Apply N Rule: no

Isotope RI (%): 1.00

MSn Logic Mode: AND


Electron Ions: both

Use MSn Info: yes


Isotope Res: 10000

Max Results: 10


Event#: 1 MS(E+) Ret. Time : 0.920 > 1.107 > 0.227 > 0.464 Scan# : 139 > 167 > 35 > 71



Measured region for 251.1155 m/z






C16 H14 N2 O [M+H]+ : Predicted region for 251.1179 m/z



| Rank | Score | Formula (M) | Ion | Meas. m/z | Pred. m/z | Df. (mDa) | Df. (ppm) | Iso | DBE |
|------|-------|--------------|--------|-----------|-----------|-----------|-----------|-------|------|
| 2 | 25.80 | C16 H14 N2 O | [M+H]+ | 251.1155 | 251.1179 | -2.4 | -9.56 | 58.11 | 11.0 |

(R)-5-(2-fluorophenyl)-3-methyl-1,3-dihydro-2H-benzo[e][1,4]diazepin-2-one(5d)



Formula Predictor Report - COMP 11_12082020 Analysis_31.lcd

Page 1 of 1

Data File: C:\LabSolutions\OldData\Viashini Rajaratnam\12082020 Analysis\COMP 11_12082020 Analysis_31.lcd

| Elmt | Val. | Min | Max | Use Adduct |
|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------------|
| H | 1 | 10 | 30 | N | 3 | 1 | 3 | Si | 4 | 0 | 0 | Br | 1 | 0 | 0 | H |
| 2H | 1 | 0 | 0 | O | 2 | 1 | 3 | S | 2 | 0 | 0 | I | 3 | 0 | 0 | K |
| C | 4 | 10 | 25 | F | 1 | 0 | 1 | Cl | 1 | 0 | 1 | | | | | |

Error Margin (ppm): 300

HC Ratio: unlimited

Max Isotopes: all

MSn Iso RI (%): 75.00

DBE Range: -100.0 - 2000.0

Apply N Rule: no

Isotope RI (%): 1.00

MSn Logic Mode: AND


Electron Ions: both

Use MSn Info: yes


Isotope Res: 10000

Max Results: 10


Event#: 1 MS(E+) Ret. Time : 0.987 > 1.507 + 0.240 > 0.420 Scan# : 149 > 227 + 37 > 63






Measured region for 269.1061 m/z



C16 H13 N2 O F [M+H]+ : Predicted region for 269.1085 m/z



| Rank | Score | Formula (M) | Ion | Meas. m/z | Pred. m/z | Dif. (mDa) | Df. (ppm) | Iso | DBE |
|------|-------|----------------|--------|-----------|-----------|------------|-----------|-------|------|
| 3 | 47.48 | C16 H13 N2 O F | [M+H]+ | 269.1061 | 269.1085 | -2.4 | -8.92 | 93.47 | 11.0 |

(S)-5-(2-fluorophenyl)-3-methyl-1,3-dihydro-2H-benzo[e][1,4]diazepin-2-one (5e)

Formula Predictor Report - COMP 10_12082020 Analysis_33.lcd

Page 1 of 1

Data File: C:\LabSolutions\Data\Vilashini Rajaratnam\12082020 Analysis\COMP 10_12082020 Analysis_33.lcd


| Elmt | Val. | Min | Max | Use Adduct |
|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------------|
| H | 1 | 10 | 30 | N | 3 | 1 | 3 | Si | 4 | 0 | 0 | Br | 1 | 0 | 0 | H |
| 2H | 1 | 0 | 0 | O | 2 | 1 | 3 | S | 2 | 0 | 0 | I | 3 | 0 | 0 | K |
| C | 4 | 10 | 25 | F | 1 | 0 | 1 | Cl | 1 | 0 | 1 | | | | | |

Error Margin (ppm): 300
 HC Ratio: unlimited
 Max Isotopes: all
 MSn Iso RI (%): 75.00


DBE Range: -100.0 - 2000.0
 Apply N Rule: no
 Isotope RI (%): 1.00
 MSn Logic Mode: AND

Electron Ions: both
 Use MSn Info: yes
 Isotope Res: 10000
 Max Results: 10


Event#: 1 MS(E+) Ret. Time : 1.533 > 1.680 + 0.067 > 1.089 Scan# : 231 > 253 + 11 > 165





Measured region for 269.1062 m/z



C16 H13 N2 O F [M+H]+ : Predicted region for 269.1085 m/z



| Rank | Score | Formula (M) | Ion | Mess. m/z | Pred. m/z | Df. (mDa) | Df. (ppm) | Iso | DBE |
|------|-------|----------------|--------|-----------|-----------|-----------|-----------|-------|------|
| 4 | 46.02 | C16 H13 N2 O F | [M+H]+ | 269.1062 | 269.1085 | -2.3 | -8.55 | 84.45 | 11.0 |

(S)-5-(2-chlorophenyl)-3-methyl-1,3-dihydro-2H-benzo[e][1,4]diazepin-2-one (5f)

Formula Predictor Report - MYMV6_12082020 Analysis_19.lcd

Page 1 of 1

Data File: C:\LabSolutions\Dat\Viashini Rajaratnam\12082020 Analysis\MYMV6_12082020 Analysis_19.lcd

| Elmt | Val. | Min | Max | Use Adduct |
|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------------|
| H | 1 | 10 | 30 | N | 3 | 1 | 3 | Si | 4 | 0 | 0 | Br | 1 | 0 | 0 | H |
| 2H | 1 | 0 | 0 | O | 2 | 1 | 3 | S | 2 | 0 | 0 | I | 3 | 0 | 0 | K |
| C | 4 | 10 | 25 | F | 1 | 0 | 0 | Cl | 1 | 0 | 1 | | | | | |

Error Margin (ppm): 300

HC Ratio: unlimited

Max Isotopes: all

MSn Iso RI (%): 75.00

DBE Range: +100.0 - 2000.0

Apply N Rule: no

Isotope RI (%): 1.00

MSn Logic Mode: AND


Electron Ions: both

Use MSn Info: yes


Isotope Res: 10000

Max Results: 10


Event#: 1 MS(E+) Ret. Time : 1.253 > 1.427 > 0.120 > 0.471 Scan# : 189 > 215 > 19 > 71



Measured region for 285.0767 m/z






C16 H13 N2 O Cl [M+H]+ : Predicted region for 285.0789 m/z



| Rank | Score | Formula (M) | Ion | Meas. m/z | Pred. m/z | Df. (mDa) | Df. (ppm) | Iso | DBE |
|------|-------|-----------------|--------|-----------|-----------|-----------|-----------|-------|------|
| 1 | 62.41 | C16 H13 N2 O Cl | [M+H]+ | 285.0767 | 285.0789 | -2.2 | -7.72 | 99.38 | 11.0 |

(S)-3-methyl-5-(pyridin-2-yl)-1,3-dihydro-2H-benzo[e][1,4]diazepin-2-one(5g)



Formula Predictor Report - COMP 18_12082020 Analysis_21.lcd

Page 1 of 1

Data File: C:\LabSolutions\OldData\Viashini Rajaratnam\12082020 Analysis\COMP 18_12082020 Analysis_21.lcd


| Elmt | Val. | Min | Max | Use Adduct |
|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------------|
| H | 1 | 10 | 30 | N | 3 | 1 | 3 | Si | 4 | 0 | 0 | Br | 1 | 0 | 0 | H |
| 2H | 1 | 0 | 0 | O | 2 | 1 | 3 | S | 2 | 0 | 0 | I | 3 | 0 | 0 | K |
| C | 4 | 10 | 25 | F | 1 | 0 | 0 | Cl | 1 | 0 | 1 | | | | | |

Error Margin (ppm): 300
HC Ratio: unlimited
Max Isotopes: all
MSn Iso RI (%): 75.00


DBE Range: +100.0 - 2000.0
Apply N Rule: no
Isotope RI (%): 1.00
MSn Logic Mode: AND

Electron Ions: both
Use MSn Info: yes
Isotope Res: 10000
Max Results: 10


Event#: 1 MS(E+) Ret. Time : 1.133 > 1.507 > 0.253 > 0.537 Scan# : 171 > 227 > 39 > 81





Measured region for 252.1113 m/z



C15 H13 N3 O [M+H]+ : Predicted region for 252.1131 m/z



| Rank | Score | Formulas (M) | Ion | Meas. m/z | Pred. m/z | Df. (mDa) | Df. (ppm) | Iso | DBE |
|------|-------|--------------|--------|-----------|-----------|-----------|-----------|-------|------|
| 1 | 43.70 | C15 H13 N3 O | [M+H]+ | 252.1113 | 252.1131 | +1.8 | +7.14 | 63.70 | 11.0 |

5-(2-fluorophenyl)-2-oxo-2,3-dihydro-1H-benzo[e][1,4]diazepin-3-yl acetate(5h)

Formula Predictor Report - COMP 13_12082020 Analysis_27.lcd

Page 1 of 1

Data File: C:\LabSolutions\Dat\Viashini Rajaratnam\12082020 Analysis\COMP 13_12082020 Analysis_27.lcd

| Elmt | Val. | Min | Max | Use Adduct |
|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------------|
| H | 1 | 10 | 30 | N | 3 | 1 | 3 | Si | 4 | 0 | 0 | Br | 1 | 0 | 0 | H |
| 2H | 1 | 0 | 0 | O | 2 | 1 | 3 | S | 2 | 0 | 0 | I | 3 | 0 | 0 | K |
| C | 4 | 10 | 25 | F | 1 | 0 | 1 | Cl | 1 | 0 | 1 | | | | | |

Error Margin (ppm): 300

HC Ratio: unlimited

Max Isotopes: all

MSn Iso RI (%): 75.00

DBE Range: -100.0 - 2000.0

Apply N Rule: no

Isotope RI (%): 1.00

MSn Logic Mode: AND


Electron Ions: both

Use MSn Info: yes


Isotope Res: 10000

Max Results: 10


Event#: 1 MS(E+) Ret. Time : 1.147 -> 1.480 -> 0.240 -> 0.600 Scan# : 173 -> 223 -> 37 -> 91






Measured region for 313.0964 m/z



C17 H13 N2 O3 F [M+H]+ : Predicted region for 313.0983 m/z



| Rank | Score | Formula (M) | Ion | Meas. m/z | Pred. m/z | Df. (mDa) | Df. (ppm) | Iso | DBE |
|------|-------|-----------------|--------|-----------|-----------|-----------|-----------|-------|------|
| 3 | 74.08 | C17 H13 N2 O3 F | [M+H]+ | 313.0964 | 313.0983 | +1.9 | +6.07 | 93.42 | 12.0 |

(2-aminophenyl)(phenyl)methanone (7a)

Formula Predictor Report - COMP 1_12082020 Analysis_51.lcd

Page 1 of 1

Data File: C:\LabSolutions\LabData\Vilashini Rajaratnam\12082020 Analysis\COMP 1_12082020 Analysis_51.lcd

| Elmt | Val. | Min | Max | Use Adduct |
|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------------|
| H | 1 | 10 | 15 | N | 3 | 1 | 3 | Si | 4 | 0 | 0 | Br | 1 | 0 | 0 | H |
| 2H | 1 | 0 | 0 | O | 2 | 1 | 3 | S | 2 | 0 | 0 | I | 3 | 0 | 0 | K |
| C | 4 | 10 | 25 | F | 1 | 0 | 1 | Cl | 1 | 0 | 0 | | | | | |

Error Margin (ppm): 300

DBE Range: -100.0 + 2000.0

Electron Ions: both

HC Ratio: unlimited

Apply N Rule: no

Use MSn Info: yes

Max Isotopes: all

Isotope RI (%): 1.00


Isotope Res: 10000

MSn Iso RI (%): 75.00


MSn Logic Mode: AND

Max Results: 10


Event#: 1 MS(E+) Ret. Time : 0.960 > 1.213 - 0.160 > 0.193 Scan# : 145 > 183 - 25 > 29



Measured region for 198.0903 m/z






C13 H11 N O [M+H]+ : Predicted region for 198.0913 m/z



| Rank | Score | Formula (M) | Ion | Meas. m/z | Pred. m/z | Dif. (mDa) | Dif. (ppm) | Iso | DBE |
|------|-------|-------------|--------|-----------|-----------|------------|------------|-------|-----|
| 1 | 53.40 | C13 H11 N O | [M+H]+ | 198.0903 | 198.0913 | +1.0 | +5.05 | 59.66 | 9.0 |

(2-aminophenyl)(2-fluorophenyl)methanone (7b)



Formula Predictor Report - COMP 3_12082020 Analysis_47.lcd

Page 1 of 1

Data File: C:\LabSolutions\Datas\ Vilashini Rajaratnam\12082020 Analysis\COMP 3_12082020 Analysis_47.lcd

| Elmt | Val. | Min | Max | Use Adduct |
|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------------|
| H | 1 | 10 | 30 | N | 3 | 1 | 3 | Si | 4 | 0 | 0 | Br | 1 | 0 | 0 | H |
| 2H | 1 | 0 | 0 | O | 2 | 1 | 3 | S | 2 | 0 | 0 | I | 3 | 0 | 0 | K |
| C | 4 | 10 | 25 | F | 1 | 0 | 1 | Cl | 1 | 0 | 0 | | | | | |

Error Margin (ppm): 300

DBE Range: -100.0 - 2000.0

Electron Ions: both

HC Ratio: unlimited

Apply N Rule: no

Use MSn Info: yes

Max Isotopes: all

Isotope RI (%): 1.00


Isotope Res: 10000

MSn Iso RI (%): 75.00


MSn Logic Mode: AND

Max Results: 10


Event#: 1 MS(E+) Ret. Time : 1.480 > 1.533 + 0.200 > 0.725 Scan# : 223 > 231 + 31 > 109



Measured region for 216.0798 m/z






C13 H10 N O F [M+H]+ : Predicted region for 216.0819 m/z



| Rank | Score | Formula (M) | Ion | Meas. m/z | Pred. m/z | Df. (mDa) | Df. (ppm) | Iso | DBE |
|------|-------|---------------|--------|-----------|-----------|-----------|-----------|-------|-----|
| 2 | 33.60 | C13 H10 N O F | [M+H]+ | 216.0798 | 216.0819 | -2.1 | -9.72 | 78.51 | 9.0 |

(2-aminophenyl)(2-chlorophenyl)methanone(7c)



Formula Predictor Report - COMP 6_12082020 Analysis_41.lcd

Page 1 of 1

Data File: C:\LabSolutions\1Data\Viashini Rajaratnam\12082020 Analysis\COMP 6_12082020 Analysis_41.lcd

| Elmt | Val. | Min | Max | Use Adduct |
|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------------|
| H | 1 | 10 | 25 | N | 3 | 1 | 3 | Si | 4 | 0 | 0 | Br | 1 | 0 | 0 | H |
| 2H | 1 | 0 | 0 | O | 2 | 1 | 3 | S | 2 | 0 | 0 | I | 3 | 0 | 0 | K |
| C | 4 | 10 | 25 | F | 1 | 0 | 1 | Cl | 1 | 0 | 1 | | | | | |

Error Margin (ppm): 300

HC Ratio: unlimited

Max Isotopes: all

MSn Iso RI (%): 75.00

DBE Range: +100.0 ~ 2000.0

Apply N Rule: no

Isotope RI (%): 1.00

MSn Logic Mode: AND


Electron Ions: both

Use MSn Info: yes


Isotope Res: 10000

Max Results: 10


Event#: 1 MS(E+) Ret. Time : 0.893 > 0.973 > 1.227 > 1.238 Scan# : 135 > 147 > 185 > 187



Measured region for 232.0506 m/z






C13 H10 N O Cl [M+H]+ : Predicted region for 232.0524 m/z



| Rank | Score | Formula (M) | Ion | Meas. m/z | Pred. m/z | Df. (mDa) | Df. (ppm) | Iso | DBE |
|------|-------|----------------|--------|-----------|-----------|-----------|-----------|-------|-----|
| 1 | 30.35 | C13 H10 N O Cl | [M+H]+ | 232.0506 | 232.0524 | -1.8 | -7.76 | 48.64 | 9.0 |

(2-aminophenyl)(pyridin-2-yl)methanone (7d)



Formula Predictor Report - COMP 7_12082020 Analysis_39.lcd

Page 1 of 1

Data File: C:\LabSolutions\Data\Vilashini Rajaratnam\12082020 Analysis\COMP 7_12082020 Analysis_39.lcd

| Elmt | Val. | Min | Max | Use Adduct |
|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------------|
| H | 1 | 10 | 30 | N | 3 | 1 | 3 | Si | 4 | 0 | 0 | Br | 1 | 0 | 0 | H |
| 2H | 1 | 0 | 0 | O | 2 | 1 | 3 | S | 2 | 0 | 0 | I | 3 | 0 | 0 | K |
| C | 4 | 10 | 25 | F | 1 | 0 | 1 | Cl | 1 | 0 | 1 | | | | | |

Error Margin (ppm): 300

DBE Range: -100.0 -> 2000.0

Electron Ions: both

HC Ratio: unlimited

Apply N Rule: no

Use MSn Info: yes

Msx Isotopes: all

Isotope RI (%): 1.00


Isotope Res: 10000

MSn Iso RI (%): 75.00


MSn Logic Mode: AND

Max Results: 10


Event#: 1 MS(E+) Ret. Time : 0.893 -> 1.640 -> 0.040 -> 0.311 Scan# : 135 -> 247 -> 7 -> 47



Measured region for 199.0844 m/z





C12 H10 N2 O [M+H]+ : Predicted region for 199.0866 m/z



| Rank | Score | Formula (M) | Ion | Meas. m/z | Pred. m/z | Df. (mDa) | Df. (ppm) | Iso | DBE |
|------|-------|--------------|--------|-----------|-----------|-----------|-----------|-------|-----|
| 1 | 35.70 | C12 H10 N2 O | [M+H]+ | 199.0844 | 199.0866 | -2.2 | -11.05 | 95.97 | 9.0 |

ethyl (R)-6-(2-fluorophenyl)-4-methyl-4H-benzo[f]imidazo[1,5-a][1,4]diazepine-3-carboxylate(9a)



Formula Predictor Report - COMP 14_12082020 Analysis_25.lcd

Page 1 of 1

Data File: C:\LabSolutions\1Data\Nilashini Rajarathnam\12082020 Analysis\COMP 14_12082020 Analysis_25.lcd


| Elmt | Val. | Min | Max | Elmt | Val. | Min | Max | Elmt | Val. | Min | Max | Use Adduct |
|------|------|-----|-----|------|------|-----|-----|------|------|-----|-----|------------|
| H | 1 | 10 | 30 | O | 2 | 1 | 3 | Cl | 1 | 0 | 0 | H |
| 2H | 1 | 0 | 0 | F | 1 | 1 | 4 | Br | 1 | 0 | 0 | Na |
| C | 4 | 10 | 25 | Si | 4 | 0 | 0 | I | 3 | 0 | 0 | K |
| N | 3 | 1 | 3 | S | 2 | 0 | 0 | | | | | NH4 |

Error Margin (ppm): 300
 HC Ratio: unlimited
 Max Isotopes: all
 MSn Iso RI (%): 75.00


DBE Range: -100.0 - 2000.0
 Apply N Rule: no
 Isotope RI (%): 1.00
 MSn Logic Mode: AND

Electron Ions: both
 Use MSn Info: yes
 Isotope Res: 10000
 Max Results: 10


Event#: 1 MS(E+) Ret. Time : 0.907 > 1.133 - 0.320 > 0.467 Scan# : 137 > 171 - 49 > 71



Measured region for 364.1437 m/z




C21 H18 N3 O2 F [M+H]+ : Predicted region for 364.1456 m/z



| Rank | Score | Formula (M) | Ion | Meas. m/z | Pred. m/z | Df. (mDa) | Df. (ppm) | Iso | DBE |
|------|-------|-----------------|--------|-----------|-----------|-----------|-----------|-------|------|
| 3 | 72.18 | C21 H18 N3 O2 F | [M+H]+ | 364.1437 | 364.1456 | +1.9 | +5.22 | 82.21 | 14.0 |


Debromination on a double bond containing substrate




PHMS debrominated the bromine group selectively without interfering with the double bond. The debromination was confirmed by running the crude products on a single quadrupole LCMS 2020.

<Spectrum>


R.Time:----(Scan#:)----
 MassPeaks:221 BasePeak:118(784464)
 Spectrum Mode:Averaged 0.093-0.313(29-95)
 BG Mode:None Polarity:Positive Segment 1 - Event 1




R.Time:----(Scan#:)----
 MassPeaks:231 BasePeak:91(194980)
 Spectrum Mode:Averaged 0.097-0.317(30-96)
 BG Mode:None Polarity:Negative Segment 1 - Event 2



Prithu 2020 - 3-87/3-526-1 - 3.lcd




If double bond is reduced, then the mass will be 120, there is no peak here.



Exact Mass: 120.09

If double bond is reduced before debromination, then the mass will be 198 and a double peak will appear but there is no such peak in the whole spectra



Exact Mass: 198.00